


Validation of BacT/ALERT® Media for Sterility Testing in a Skan PSI Isolator after Three H₂O₂ Decontamination Cycles

Mary Beth Anheuser¹, Daniel Eichenberger³, Diana Tehsmann³, Carolyn Beres¹, Michael Ullery¹, Scott Jeffrey², Jessica Cook¹, Gregory Devulder¹
¹bioMérieux, Inc. Hazelwood, MO / ²bioMérieux, Inc. Durham, NC / ³SKAN AG, Switzerland



INTRODUCTION

Pharmaceutical and biotechnology manufacturers using the BacT/ALERT® Dual-T (BTA) system for sterility testing of aseptically produced products conduct the test under aseptic conditions to mitigate the risk of false positive results from environmental contamination. The Skan PSI Isolator, using H₂O₂ for decontamination, provides an aseptic environment where samples can be introduced into the BTA media. Typically, enough media are decontaminated to perform a scheduled test. Unused media are removed from the isolator and discarded unless performance is validated after multiple decontamination cycles. Additionally, interruptions to the decontamination cycle such as an electrical event can result in additional H₂O₂ exposure to the media and potentially to laboratory personnel.

A study was performed to validate BTA media for sterility testing after exposure to three isolator decontamination cycles.

MATERIALS AND METHODS

Decontamination Cycle:

BTA media (iAST, iNST, iFA, iPF, iFN, iLYM) and water-filled bottles (1 ml purified) were exposed to three consecutive decontamination cycles in a Skan Pharmaceutical Isolator (PSI), type SIS-700, equipped with a vaporizer.

Set Values of the cycle parameters:

Parameter	Set Value
q 1 Quantity of H ₂ O ₂ vaporized during main dosage (g)	70
q 2 Quantity of H ₂ O ₂ vaporized during redosage (g)	7
x Number of redosages	2
t 1 Holding time between redosage (sec)	180
t 2 Aeration time (sec)	10800
Relative humidity of start conditioning phase (%)	< 20
Temperature of start conditioning phase (°C)	< 40
Concentration of the H ₂ O ₂ solution (%)	49.5

Top caps on the iNST, iFA, iPF and iFN bottles were removed prior to decontamination. Two load configurations were tested: 209 bottles and 409 bottles over six shelves of two racks. Three biological indicators (*Geobacillus stearothermophilus*, NAS-152, Apex Laboratories, Sanford, NC) were placed during the third cycle to validate that decontamination was achieved for each load pattern tested.

Chemical-Physical Testing:

Residual H₂O₂ was measured using the Amplex® UltraRed Assay (Life Technologies™). Nine decontaminated bottles from each load configuration and an equal number of controls were rinsed with 10 ml of purified water and tested for residual H₂O₂ (average H₂O₂ (mg) per bottle surface). Penetration of H₂O₂ was determined by direct testing of the water recovered from each of nine bottles per load configuration and an equal number of controls (average H₂O₂ ng / ml). The closures on 100 decontaminated and 100 control bottles were tested using a 1600 Type M durometer (Rex Gauge Co., Buffalo Grove, IL) to determine if exposure to H₂O₂ caused a change in the hardness of the stopper. All bottles were inspected for evidence that exposure to H₂O₂ had no impact on container integrity.

Growth Promotion Testing:

The growth promoting properties of control and decontaminated media were tested using BioBall® (BTF, Sydney, Australia), 25 cfu / bottle, five replicates of each microorganism / medium / incubation temperature combination listed in Table 1: *Staphylococcus aureus* NCTC10788 (SA), *Escherichia coli* NCTC12923 (EC), *Clostridium sporogenes* NCTC12935 (CS), *Bacillus subtilis* NCTC10400 (BS), *Pseudomonas aeruginosa* NCTC12924 (PA), *Candida albicans* NCPF3179 (CA), and *Aspergillus brasiliensis* NCPF2275 (AB).

RESULTS & DISCUSSION

GROWTH PROMOTION TESTING (GPT):

The chemical and physical test results were correlated to the GPT results. If analysis of Time to Detection (TTD) confirmed the media remained unaffected after decontamination, any differences in the chemical and physical tests were considered insignificant. A difference in mean TTD between the test and control bottles was considered statistically significant if p-value = <0.05. A difference between means with a positive value is interpreted as a decrease in TTD for the test bottles and a negative value is interpreted as an increase in TTD for the test bottles.

Table 1: Growth Promotion Test Results

Medium	T°C	Mean TTD, Hours				Difference, Hours	P-Value	
		Cntrl	STD	Test	STD			
SA	32.5°C	i AST	23.4	0.6	24.9	1.1	-1.5	0.025*
		i FA	21.8	1	21.9	1.2	-0.1	0.869
		i PF	20.8	0.7	21.4	0.4	-0.6	0.154
		i FN	23	0.8	23.5	0.6	-0.5	0.320
EC	32.5°C	i NST	19.5	0.8	19.9	0.9	-0.3	0.541
		i AST	16.7	0.4	16.5	0.5	0.2	0.514
		i FA	17.7	0.3	16.6	0.6	1.1	0.008*
		i PF	23.7	0.6	22.3	0.8	1.4	0.014*
CS	32.5°C	i FN	15.4	0.3	15	0.4	0.4	0.086
		i NST	14.6	0.1	14.4	0.3	0.2	0.252
		i FN	23	1.2	23.8	0.7	-0.8	0.232
		i NST	22.1	1.2	21.7	0.7	0.4	0.565
BS	32.5°C	i AST	18.4	0.3	18.1	0.7	0.2	0.538
		i FA	19.4	0.4	19.3	0.3	0.1	0.613
		i PF	18.6	0.4	18.9	0.7	-0.3	0.408
		i AST	41.5	2.5	43.2	2.4	-1.7	0.317
PA	22.5°C	i FA	44.9	2.8	44.9	0.8	0	0.977
		i PF	45.1	2.2	45.8	2.8	-0.8	0.656
		i AST	22.8	0.7	22.5	0.5	0.3	0.482
		i FA	22.2	0.5	22.5	0.4	-0.3	0.357
CA	32.5°C	i PF	22.1	0.4	22.2	0.6	-0.1	0.665
		i AST	45	1.1	43.6	0.4	1.4	0.028*
		i FA	44.9	3.1	43.5	0.8	1.3	0.404
		i PF	43.3	1	43.8	1	-0.5	0.471
AB	32.5°C	i AST	28.6	0.7	28.5	0.8	0	0.933
		i FA	44.2	1.6	42.3	2.5	2	0.182
		i PF	56.2	6	58.6	6.9	-2.4	0.575
		i LYM	25.2	0.5	25.5	1	-0.3	0.555
AB	22.5°C	i AST	52.8	0	49	2.1	3.8	0.016*
		i FA	53.8	2.1	50.9	1.1	2.9	0.028*
		i LYM	46.5	1	45.1	0.8	1.4	0.004*
		i PF	56.6	1.3	52.3	2	4.3	0.034*
AB	32.5°C	i AST	58.6	2.7	58.1	3.1	0.5	0.803
		i FA	58.6	10.1	55.2	5.1	3.4	0.526
		i PF	51.8	2.1	51.8	1.3	0	1.000
		i LYM	32.5	2.1	33.4	2.1	-0.9	0.540
AB	22.5°C	i AST	106.6	3.1	104.6	2.5	1.9	0.641
		i FA	102.2	9.1	102.7	7.3	-0.5	0.929
		i PF	96.5	4.3	105.1	10	-8.6	0.113
		i LYM	59.5	2	60.5	3.6	-1	0.614

*Indicates statistically significant difference in TTD.

Eight microorganism / medium / incubation temperature combinations had p-values <0.05.

Although statistically significant, the results were considered practically insignificant since seven of eight results for the decontaminated bottles showed improvements in TTD. The GPT results indicate that the decontaminated bottle H₂O₂ penetration levels were also insignificant.

CHEMICAL-PHYSICAL TESTING:

Negative Control: Biological indicators from each load were negative for growth after 7 d incubation at 55-60 C. This validated that the cycle parameters were effective for Decontamination of the load configurations tested.

Container – Stopper: No visible effects to the container were observed and label and barcode readability were maintained. Residual amounts of H₂O₂ found on the surface of the decontaminated bottles were measured to provide information on potential user exposure (Table 2). Durometer hardness measurements show that three H₂O₂ decontamination cycles do not cause the bottle stopper to harden and impact its ability to seal (Table 4).

Residues: Mean levels of H₂O₂ inside the decontaminated bottles were higher than the controls, however, the levels were still very low – 17.15 and 6.16 ng / ml (ppb) (Table 3).

Table 2: Surface Residual levels of H₂O₂

Surface Residual	Load 1		Load 2	
	Mean, mg / bottle	Std Dev	Mean, mg / bottle	Std Dev
Test	0.9152	0.2379	1.3402	0.2547
Control	0.0001	0.0001*	0.0002	0.0002
Total	0.9151		1.34	

*4 of 9 control bottles were below the detection limit of the test (3.9 ng H₂O₂/ml)

Table 3: Bottle Penetration levels of H₂O₂

Penetration	Load 1		Load 2	
	Mean, ng / mL	Std Dev	Mean, ng / mL	Std Dev
Test	47.8	12.07	17.08	11.17*
Control	30.65	3.97	10.92	5.62*
Total	17.15		6.16	

*1 of 9 bottles was below the detection limit of the test (3.9 ng H₂O₂/ml)

Table 4: Durometer Hardness

Stopper	Durometer Hardness Measurement, Shore M			
	Mean	Std Dev	Minimum	Maximum
Test	53.08	2.74	48	58
Control	53.87	2.34	49	58

CONCLUSION

This study verified the integrity of the BTA media container and validated that the BTA media remained unaffected after exposure to three typical H₂O₂ decontamination cycles. BacT/ALERT bottles (iAST, iNST, iFA, iFN, iPF and iLYM) are suitable for use in an Isolator.

REFERENCES

- Guidance for Industry: Sterile Drug Products Produced by Aseptic Processing-Current Good Manufacturing Process (September 2004).
- United States Pharmacopeia <1208> Sterility Testing-Validation of Isolator Systems.
- United States Pharmacopeia <71> Sterility Tests.
- ISO 1138-1 Sterilization of Healthcare Products – Biological Indicators, General.
- PIC's Recommendation "Isolators used for Aseptic Processing and Sterility Testing" PI 014-2, June 2004.
- European Pharmacopoeia 5.1.9. Guidelines For Using The Test For Sterility.
- Sigwarth, V. and Moirandat, C., "Development and Quantification of H₂O₂ Decontamination Cycles", PDA Journal of Pharmaceutical Science and Technology, July/August (2000): 286-304.
- Beers, R.F., Sizer, I.W. (1952): A spectrophotometric method for measuring the breakdown of hydrogen peroxide by catalase (J. Biol. Chem. 195: p. 133-140).