



Superior design
for efficient and safe work

BIOWIZARD *Silver Line*



BIOWIZARD
KOJAIR®

Microbiological
Safety Cabinet
class II EN 12469 certified

Kojair quality
for reasonable price



BLOWIZARD *Silver Line*

Microbiological Safety Cabinet class II certified according to EN 12469

Standard features

- UV-light outlet
- UV-timer
- Electric front window
- Large side windows
- Double electric socket
- Timer for electric socket
- Switch for electric socket
- Switch for gas valve
- Stand-by switch



Easy installation

- The cabinet will easily fit through common, 80 cm wide door openings
- The height of the working level can be chosen on site: 700, 800 or 900 mm
- The cabinet has connections for fixed or thimble type exhaust ducting (optional)

Easy service

- All maintenance can be performed from the front
- Easy performance of annual testing procedures

Alarms

- Acoustic and visual alarms
in any failure of air flows
in a false front window position





Kojair ergonomics ensure safe and practical working conditions.



Kojair ergonomics

- Silent and comfortable, sound level < 55 dBa
- Efficient light tubes ensure shadowless lighting over the whole working area
- Completely closing electrical sliding front window
- Adjustable base stand



Easy to clean

- All surfaces are smooth without any screws or crevices
- The polished stainless steel three-piece table top is easy to remove for cleaning
- The bottom sink is stainless steel with round corners
- The front window is easily cleaned both inside and outside
- Side and front windows are made of laminated glass

BIOWIZARD
KOJAIR®



Dimensions and Options

Biowizard Silver Line	100	130	170	200
Total width (mm) B	975	1280	1625	1930
Working space (mm)				
Width	885	1190	1535	1840
Depth	600	600	600	600
Height	720	720	720	720
Air volumes (m ³ /h)				
Recirculation	740	995	1320	1580
Exhaust	255	345	445	530
Total	995	1340	1670	2010
• Exhaust fitting (ø, mm)	160	160	200	200
• Exhaust fitting location A	244	288	355	733
Lighting (W)	2x40	2x40	2x40	2x40
Weight (incl. base, kg)	180	240	290	330
Power (max. W)	880	1150	1520	1540
Transport height with base 1990 mm, without base 1440 mm				



• = option

Options

- Adjustable base stand, electrical
- Gas tap with magnetic valve
- Exhaust duct connector, fixed
- Extra electric socket
- Vacuum tap
- Exhaust duct connector, thimble type
- Socket in the backchannel
- Compressed air tap
- Antiflow back valve
- Table top in one piece
- UV light with magnetics
- Potential free signals

Leading expertise in **clean air technology**



Kojair Tech Oy
Teollisuustie 3
FI-35700 Vilppula

Tel. +358 3 471 7000
Fax +358 3 471 7077

sales@kojair.com
www.kojair.com